EN Instructions for use: Purpose: Natural rubber latex condoms are intended for single use. They serve to help prevent sexually transmitted infections and pregnancy when used as directed during vaginal intercourse. Natural rubber latex condoms are transparent, with a reservoir tip and lubricant. Intended User Population: This product is intended for use by adult males as a contraceptive and prophylactic measure. Instructions: Hold the reservoir at the tip of the condom between your thumb and forefinger to prevent any air pockets from forming. Roll the condom down over the erect penis. The condom must be put on before the first contact between the penis and your partner's body to help prevent sexually transmitted infections and pregnancy. Withdraw the penis soon after ejaculation, holding the condom at the base. Contraindications: Do not use this product if you have a known allergy to natural rubber latex. Undesirable side effects: Allergic reactions may occur in sensitive individuals, including itching, redness, rash, or irritation. In rare cases, severe allergic reactions such as anaphylactic shock may occur. Condom breakage or leakage may rarely occur, which may increase the risk of unintended pregnancy or STIs. Proper storage and usage can minimize these risks. Warnings: If you feel the condom slipping off or becoming too tight on the penis, stop and check the fit, as this could lead to tearing. You may need to choose a different condom size. If using additional lubricant, it is essential to use only those specifically designed for use with condoms. Avoid using oil-based lubricants such as petroleum jelly, baby oil, body lotions, massage oils, butter, margarine, etc., as these can compromise the condom's integrity. If you are using any medication that may come into contact with the condom, consult a doctor or pharmacist about its compatibility with the condom. Use the condom only once, as reusing it increases the risk of failure or infection. If a condom leaks or breaks during use, seek medical help as soon as possible, and at least within 72 hours. Do not use the product if the rubber material is sticky or brittle, has cracks, or is visibly damaged. **No contraceptive method can guarantee 100% protection.** Precautions: Handle the condom carefully when removing it from the packaging to avoid damage from fingernails, jewelry, etc. If the individual wrapper (foil) is visibly damaged, discard that condom and use a new one from an undamaged package. Dispose of the used condom in the trash, not in the toilet. Additional information: Please note the batch number behind the batch code symbol and the expiry date behind the "use by" symbol on the packaging. Please store the condoms in their original packaging, in a cool, dry place, and not outside the temperature range indicated by the warning symbol or in unsuitable locations (e.g., car, wallet, etc.). Short term exposure to temperature in between –15°C to +50°C is permitted during transportation. These condoms comply with the international standards ISO 4074 / EN ISO 4074 and are individually electrically tested. Incident Reporting: Any serious incident related to the product must be reported to the manufacturer and the competent authority in the Member State where the user resides.

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